New Brunswick™ TCA-3 Temperature Monitoring System

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Copyright

Copyright © 2014 Eppendorf AG, Germany. No part of this publication may be reproduced without the prior permission of the copyright owner.

The company reserves the right to change information in this document without notice. Updates to information in this document reflect our commitment to continuing product development and improvement.

Trademarks

Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany.

New Brunswick™ is a trademark of Eppendorf AG, Germany.

VeriSign Secured® is a registered trademark of VeriSign, Inc., USA.

Microsoft® and Excel® are registered trademarks of Microsoft Corporation in the United States and other countries.

Trademarks are not marked in all cases with ™ or ® in this manual.

Eppendorf has attempted to identify the ownership of all trademarks from public records. Any omissions or errors are unintentional.

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4.4.2 In den Account einloggen

The screenshot shows the login interface for New Brunswick. On the left, there is a logo and the text 'New Brunswick an eppendorf company'. Below this, there are four input fields: 'Gateway Mac:', 'Login Name:', 'Password:', and a checkbox labeled 'Remember Gateway Mac'. A 'Login' button is positioned below the checkbox. At the bottom left, there is a link: 'Click [HERE](#) to try the demo site.' On the right side of the page, the text 'Monitoring made easy.' is displayed.

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
The screenshot shows a web page with a blue header bar containing the text "Welcome to your On". Below the header, there is a section titled "Account Setup Required" in blue text. The main text reads: "To properly setup your account we need you to complete some basic". This is followed by the sentence "In the following steps we will ask you to enter information for:" and a numbered list: "1. Primary Contact Info", "2. User Accounts Logins", and "3. The Default Alarm Contact(s)". A bold note states: "NOTE: This information is required so that we can contact you and securely configured." At the bottom right of the content area, there is a button with a dotted border and the text "Next >>".

Welcome to your On

Account Setup Required

To properly setup your account we need you to complete some basic

In the following steps we will ask you to enter information for:

1. Primary Contact Info
2. User Accounts Logins
3. The Default Alarm Contact(s)

NOTE: This information is required so that we can contact you and securely configured.

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Account Se

Primary Contact

From time to time we may need to notify you of system updates, schedules, and other important information. Who is the primary contact for this monitoring system account?

Who is the primary contact for this monitoring system account?

Full Name*:	<input type="text"/>
Company*:	<input type="text"/>
Address Line1*:	<input type="text"/> <small>Street address, P.O. box, company name, c/o</small>
Address Line2:	<input type="text"/> <small>Apartment, suite, unit, building, floor, etc.</small>
City*:	<input type="text"/>
State/Province/Region*:	<input type="text"/>
ZIP/Postal Code*:	<input type="text"/>
Country*:	<input type="text" value="United States"/>
Timezone*:	<input type="text" value="(GMT-05:00) Eastern Standard Time (DST)"/>
Phone Number*:	<input type="text"/>
Email Address*:	<input type="text"/>
Repeat Email Address*:	<input type="text"/>

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Account Se

User Accounts

Your Admin and Reader accounts have been setup with default passw
must change these passwords, if you haven't already.

<p style="margin: 0;">Administrator</p> <p style="margin: 5px 0;">Login Name*: <input style="width: 80%;" type="text" value="Admin"/></p> <p style="margin: 5px 0;">Password: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">Repeat Password: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">First Name: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">Last Name: <input style="width: 80%;" type="text"/></p>	<p style="margin: 0;">Read-Only Account</p> <p style="margin: 5px 0;">Login Name*: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">Password: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">Repeat Password: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">First Name: <input style="width: 80%;" type="text"/></p> <p style="margin: 5px 0;">Last Name: <input style="width: 80%;" type="text"/></p>
--	--

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Account Set

Default Alarm Contact

To ensure system reliability, at least one of the following lists must

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Account Set

Default Alarm Contact

To ensure system reliability, at least one of the following lists must be populated.

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
The screenshot shows the New Brunswick TCA-3 Temperature Monitoring System interface. It features a header with the company logo and navigation links (SETTINGS | SUPPORT | LOGOUT). Below the header, there is a 'Welcome' message and a 'Time to Next Data Storage' indicator. The main content is divided into two primary sections: 'TRIGGERED ALARMS' and 'LATEST MEASUREMENTS'. The 'TRIGGERED ALARMS' section contains a table with columns for Alarm Name, POD, Measurement, Value, Alarm Time, and Chart. The 'LATEST MEASUREMENTS' section shows a table with columns for Measurement Name, Value, and status icons. Red circles and lines highlight specific areas: 1. The navigation links (SETTINGS | SUPPORT | LOGOUT). 2. The 'SEND TO HISTORY' button. 3. The 'LATEST MEASUREMENTS' table. 4. The 'VIEW ALARM HISTORY' button.

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New Brunswick
an eppendorf company

SETTINGS | SUPPORT | LOGOUT

Welcome Time to Next Data Storage: 00:00:43

TRIGGERED ALARMS

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
<input type="checkbox"/> Low supply voltage	GW-0014a002017d	Supply Voltage (Battery)	2.28	01:32:00 PM 01/25/10	
<input type="checkbox"/> Power or Internet Outage	N/A	N/A		01:34:00 PM 01/25/10	

Select: All, None, Inactive ■ Active Alarm ■ Inactive Alarm

[SEND TO HISTORY](#) [VIEW ALARM HISTORY](#)

LATEST MEASUREMENTS

GW-0014a002017d	02/04/10 09:51:00 AM	RTD #1	6.58 C	RTD #2	-20.22 C

Graphs, Controls and Diagnostics | Alarms | Preferences

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


The screenshot displays the New Brunswick web interface. At the top, there is a navigation bar with "SETTINGS | SUPPORT | LOGOUT" and a "Welcome" message. Below this, a "TRIGGERED ALARMS" table is shown with the following data:

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
<input type="checkbox"/> Low supply voltage	GW-0014a002017d	Supply Voltage (Battery)	2.28	01:32:00 PM 01/25/10	
<input type="checkbox"/> Power or Internet Outage	N/A	N/A		01:34:00 PM 01/25/10	

Below the table, there are buttons for "SEND TO HISTORY" and "VIEW ALARM HISTORY". A secondary window titled "GATEWAY COMMUNICATION ALARM" is open, showing configuration options for alarm "GW-0014a002017d". The "General" tab is active, showing:

- Alarm Enabled
- Alarm Description: Name: Power or Internet Outage
- Notify List: Default List (selected), Custom List

Buttons for "OK", "Cancel", and "Apply" are at the bottom of the configuration window.

4.10 Sensoreinstellungen

The screenshot displays the New Brunswick TCA-3 Temperature Monitoring System interface. At the top, there is a navigation bar with 'SETTINGS | SUPPORT | LOGOUT' and a 'Welcome' message. Below this, the 'TRIGGERED ALARMS' section shows a table with columns: ALARM NAME, POD, MEASUREMENT, VALUE, ALARM TIME, and CHART. A single alarm is listed: 'low' for POD 'GW- 0014a002016a', measurement 'RTD #1', value '-245.69', and alarm time '08:58:00 AM 01/05/12'. Below the table are buttons for 'SEND TO HISTORY' and 'VIEW ALARM HISTORY'. The 'LATEST MEASUREMENTS' section shows a table with columns: Measurement Name, Value, and three action icons (Chart, Alarm, Preference). The table lists: RTD #1 (-245.69 C), RTD #2 (-237.56 C), Supply Voltage (Battery) (4.05 Volts), Supply Voltage (External) (4.73 Volts), and Thermocouple (1,177.49 C). The 'Preference' icon for RTD #1 is circled in red. Navigation links 'Configure Data Logs and Downloads', 'Alarms', and 'Preferences' are visible above the table.

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
The screenshot shows the 'SENSOR PREFERENCES' dialog box for 'RTD #1'. The dialog has three tabs: 'General', 'Info', and 'Advanced'. The 'General' tab is active. It contains the following fields:

- Name: RTD #1
- Units: C (dropdown menu)
- Notes: (empty text area)

At the bottom of the dialog are three buttons: 'OK', 'Cancel', and 'Apply'.

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
4.11 Sensoralarme

New Brunswick
an eppendorf company

SETTINGS | SUPPORT | LOGOUT

Welcome Time to Next Data Storage: 00:01:17

TRIGGERED ALARMS

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
There are no triggered alarms					

Select: All, None, Inactive ■ Active Alarm ■ Inactive Alarm

SEND TO HISTORY VIEW ALARM HISTORY

LATEST MEASUREMENTS

Measurement	Value	Chart	Alarm	Preference
GW- 0014a002016a	01/05/12 09:19:00 AM			
RTD #1	-245.76 C			<input checked="" type="checkbox"/>
RTD #2	-237.61 C			<input checked="" type="checkbox"/>
Supply Voltage (Battery)	4.05 Volts			<input checked="" type="checkbox"/>
Supply Voltage (External)	4.73 Volts			<input checked="" type="checkbox"/>
Thermocouple	1,177.51 C			<input checked="" type="checkbox"/>

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The screenshot shows a software window titled "ALARM DEFINITION" for sensor "GW- 0014a002017d - RTD #1". It has two tabs: "General" and "Trigger". In the "General" tab, there is a checked checkbox for "Alarm Enabled". Below that is the "Alarm Description" section, which includes a "Name:" field containing the text "Fridge" and a blue "Add Notes" link. Underneath is the "Notify List:" dropdown menu, which is currently set to "Default List". A small text box below the dropdown explains: "This will notify everyone on the 'Default Notify List'. The list is located in the 'Settings' menu on the main page." At the bottom of the window are three buttons: "OK", "Cancel", and "Apply".

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New Brunswick ALARM DEFINITION

GW- 0014a002016a - RTD #1

General Trigger Email List Phone List

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The screenshot shows a software window titled "ALARM DEFINITION" for a device labeled "GW-0014a002017d - RTD #1". It has two tabs: "General" and "Trigger". The "Trigger" tab is selected and contains the following settings:

- Alarm Limits:** A section with a "[Trigger Help]" link and a "* = required" note. It includes two "Sensor value*" fields. The first has a dropdown arrow pointing left and the value "20". The second has a dropdown arrow pointing right and the value "30".
- And/Or:** A dropdown menu currently set to "OR".
- Advanced:** A section containing:
 - Trigger Mode:** A dropdown menu set to "Auto-Rearm".
 - Filter:** A field with the number "3", followed by the text "occurrences in", a field with the number "5", and the text "consecutive samples".

At the bottom of the dialog are three buttons: "OK", "Cancel", and "Apply".

4.12 Sensorkurve

The screenshot displays the New Brunswick TCA-3 Temperature Monitoring System interface. At the top, there is a navigation bar with 'SETTINGS | SUPPORT | LOGOUT' and a 'Welcome' message. Below the navigation bar, the 'TRIGGERED ALARMS' section shows a table with columns: ALARM NAME, POD, MEASUREMENT, VALUE, ALARM TIME, and CHART. The table is currently empty, with a message 'There are no triggered alarms'. Below the table, there are filters for 'Active Alarm' and 'Inactive Alarm', and buttons for 'SEND TO HISTORY' and 'VIEW ALARM HISTORY'. The 'LATEST MEASUREMENTS' section shows a table with columns: Sensor Name, Value, and three icons (Sensorkurve, Alarm, Preference). The table contains the following data:

Sensor Name	Value	Sensorkurve	Alarm	Preference
GW- 0014a002016a	01/05/12 09:19:00 AM	[Icon]	[Icon]	[Icon]
RTD #1	-245.76 C	[Icon]	[Icon]	[Icon]
RTD #2	-237.61 C	[Icon]	[Icon]	[Icon]
Supply Voltage (Battery)	4.05 Volts	[Icon]	[Icon]	[Icon]
Supply Voltage (External)	4.73 Volts	[Icon]	[Icon]	[Icon]
Thermocouple	1,177.51 C	[Icon]	[Icon]	[Icon]

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7 Zertifikate

Declaration of Conformity

The product named below fulfills the requirements of directives and standards listed. In the case of unauthorized modifications to the product or an unintended use this declaration becomes invalid.

Product name:

TCA-3
including accessories

Product type:

Remote logging temperature monitor

Relevant directives / standards:

2004/108/EC EN 61326-1
2006/95/EC EN 61010-1
2011/65/EU

Management Board

Portfolio Management

Date: January 06, 2014

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

MI288-9999-00

Eppendorf® and the Eppendorf Logo are registered trademarks of Eppendorf AG, Hamburg/Germany. All rights reserved incl. graphics and pictures. Copyright 2013 © by Eppendorf AG.

Evaluate your manual

Give us your feedback.
www.eppendorf.com/manualfeedback